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 K-Jump Health Co., Ltd.
No.56, Wu Kung 5th Rd.
New Taipei Industrial Park
New Taipei City, 24890 Taiwan

 Polygreen Germany GmbH
Ruhlsdorfer Str. 95
14532 Stahnsdorf, Germany

Importeur
SCALA Electronic GmbH
Ruhlsdorfer Str. 95
14532 Stahnsdorf, Germany

GB Classic waterproof clinical thermometer SC 1493 with flexible measuring point

CE 0197

Special features

The SC 1493 is a waterproof, digital fever thermometer with a flexible measuring tip. It displays the highest temperature within 10 seconds. It has a memory function to store the last measured value. A fever alarm beeps when the temperature measured exceeds 37.5°C.


Safety information

- The SC 1493 is suitable only for measuring human body temperatures.
- The thermometer contains small parts (batteries, etc) which could be swallowed by small children. Please keep the thermometer and its batteries out of the reach of small children. If the battery or other parts are swallowed you should seek medical help immediately.
- The enclosed batteries may not be recharged or reactivated by other methods, may not be dismantled, thrown into fires or short-circuited.
- Protect the SC 1493 from direct sunshine and temperatures above 50°C.
- Never drop the SC 1493, do not bend or open it (with the exception of the battery compartment for battery replacement).


Microbial contamination

Because a hundred percent disinfection is not possible under conditions at home, it is urgently recommended that every patient should use their own fever thermometer.

The device

The LC display

Operating manual

Press the 0/I button to switch on. A short signal tone shows that the device is in the "switched on" mode. A segment test takes place simultaneously (checks all display segments). After this, if the temperature measured by the temperature sensor is less than 32°C, the display shows an "L" and "°C". After around one second, the "°C" starts to blink. The „L“ stands for low (less than 32°C). The current temperature is continuously displayed during measurement, and the "°C" blinks during the entire measurement time. The measurement is completed after around 10 seconds and "°C" stops blinking. After the measurement is completed, the temperature value measured remains in the display until the SC 1493 is switched off. To turn the device off after use, press the 0/I button briefly. The signal tone will again be sounded. If the device is not switched off, a signal tone is sounded after around 3 minutes and the device then switches off automatically.

Fever alarm

If the temperature measured is 37.5°C or higher, the fever alarm beeps (short double beeps for about 4 seconds) and the LC display flashes.

Display of last measured value

If you keep the 0/I button pressed for longer than three seconds, the last measured value will be shown in the display. As soon as you release the 0/I button, the thermometer returns to measurement mode.


Measured value deviation

The measured temperatures will always deviate, depending on the measurement location. This is normal, and caused by deviations in temperature distribution on the human body. You could, for example, measure a rectal temperature of 36.5°C, an oral temperature of 36.0°C and an underarm value of 35.5°C.

Measurement locations

In the mouth (oral)

Small temperature variations in the mouth between measurements are normal. These are caused by varying temperature distribution within the mouth. Try to place the thermometer in one of the two heat pockets under the tongue, to the left and right of the tongue root. The measurement sensor must have good contact with the tissue. The subject should close the mouth and breathe normally through the nose so that the measured result is not affected by breathed air. Measuring period: approx. 10 seconds.


In the anus (rectal)

This is the most accurate method. It is particularly suitable for infants and small children. Carefully guide the tip of the thermometer approximately 1 cm into the anus. Measuring period: approx. 10 seconds.

Under the arm (axillary)

From a medical viewpoint this type of measurement leads to inaccurate values. It is therefore unsuitable for carrying out exact measurements.

Changing the battery

As soon as the  symbol appears in the display, the battery has run out and requires replacement. Always replace the discharged battery with one of the same type (LR 41 or AG3). To do this, proceed as follows: Pull the cover off the battery compartment. The battery can now be accessed.

Remove the used battery with the help of a pointed object (pincers, small screwdriver etc.) and put the new battery in. Please make sure it is the right way round.

Measurement checks

Checking the accuracy of the SC 1493 is not necessary in most cases.

Cleaning and disinfection

Please only use an alcohol pad. Alcohol pads can be purchased from any dispensing chemist. We recommend the use of Mucofit®-F2000 or Helipur®Hplus to disinfect the thermometer when it is being used in medical facilities. Both are effective against all hospital-typical pathogens and can be purchased from any dispensing chemist.


Guarantee

We give a guarantee period of 2 years from the date of purchase. We guarantee perfect design of this article. We will rectify all material or manufacturing errors during the guarantee period. If, contrary to expectations, defects are found, please send the device back to the importer. Please enclose the proof of purchase and a description of the defect so that the problem can be rectified as soon as possible. The guarantee excludes damages caused by improper handling, related to the battery or to the packaging.



Disposal note

The disposal of batteries is regulated. Spent batteries are collected at collection points provided by battery dealers and the communities. Dispose of the thermometer monitor as prescribed. i.e., at a collection point for electrical scrap. Dispose of the thermometer monitor as prescribed. i.e., at a collection point for electrical scrap.

Specifications

Type: maximum thermometer
Measuring range: 32,0°C to 43,9°C
Temperature below 32,0°C: display "L" for low.
Temperature above 43,9°C: display "H" for high.
0,1°C between 35,5°C and 42,0°C (0,2°C other range) at a surrounding temperature of 18°C to 28°C.
Measurement accuracy: liquid-crystal (LCD display).
Display: 1 no.: LR 41 or AG3.
Battery type: The symbol  shows when the battery is discharged.
Storage temperature: 0°C to 50°C.

Explanation of symbols

	Manufacturer
	CE-sign of a notified body (TÜV Rheinland)
	EU representative
	Product type BF
	Disposal notes
	Product period (month/year)
	Read the instruction manual

 K-Jump Health Co., Ltd.
No.56, Wu Kung 5th Rd.
New Taipei Industrial Park
New Taipei City, 24890 Taiwan

 Polygreen Germany GmbH
Ruhlsdorfer Str. 95
14532 Stahnsdorf, Germany

Importer
SCALA Electronic GmbH
Ruhlsdorfer Str. 95
14532 Stahnsdorf, Germany