

Elektrischer Läusekamm Electronic Lice Comb

SC 04

Das Gerät tötet nur erwachsene Läuse. Eier (Nissen) und Läuse im Larvenstadium werden nicht erreicht. Deshalb wird empfohlen, alle 2 Tage die Behandlung zu wiederholen. So ist sichergestellt, dass nach 3 bis 4 Wochen alle Läuse, die sich bei Beginn der Behandlung noch im Ei- oder Larvenstadium befanden, getötet wurden. Das Gerät bringt eine für den Menschen ungefährliche Spannung auf den Metallkamm.

The Electronic Lice Comb SC 04 works by transmitting a small electrical charge to the metal comb. When the metal comb is combed through dry hair, this electric charge will kill adult lice. The Electronic Lice Comb SC 04 does not kill the eggs and some may remain attached to the hair follicle. It is recommended that combing treatments are repeated every 2-3 days.





mit Licht /
with light

Kopflaus /
Head louse



(10-fach vergrößert /
enlarged 10 times)

 K-jump Health Co., Ltd.
No. 56, Wu Kung 5th Rd.
New Taipei Industrial Park
New Taipei City, 24890 Taiwan

 Polygreen Germany GmbH
Ruhlsdorfer Straße 95
14532 Stahnsdorf, Germany

Importeur/Importer
SCALA Electronic GmbH
Ruhlsdorfer Straße 95
14532 Stahnsdorf, Germany
www.scala-electronic.de


D
Bedienungsanleitung


GB
Instruction Manual


F
Mode d'emploi


I
Istruzioni per l'uso


E
Instrucciones para el uso


NL
Gebruiksaanwijzing


CZ
ávod k použití


SK
Návod na Používanie

SC 04

Elektrischer Läusekamm Electronic Lice Comb

SC 04

scala

scala



CE

TECHNISCHE DATEN
Betriebsbedingungen:
10°C - 40°C,
25%-80% rel. Luftfeuchte
Abmessung:
90 x 60 x 27 mm
Gewicht:
42 g (ohne Batterie)
Stromversorgung:
1 x AAA Alkalie-Batterie (LR03)

TECHNICAL SPECIFICATIONS
Operating environment:
10°C - 40°C, with relative
humidity: 25%-80%
Dimensions:
90 x 60 x 27 mm
Weight:
42 g (without battery)
Battery:
1 x AAA alkaline battery (LR03)

- Batterie inklusive / battery included
- batteriebetrieben / battery powered
- beleuchteter Kamm / LED Spotlight

SC 04



REV.20211013