

DE
Fieberthermometer

GB
Thermometer

FR
Thermomètre

NL
Thermometer

IT
Termometro

TR
Termometre

SC 53 FH



Ohr- und Stirn-Thermometer SC 53 FH

- ▶ Infrarot
- ▶ Messung im Ohr oder an der Stirn



0197

SC 53 FH

AVITA Corporation
9F, No. 78, Sec. 1, Kwang-Fu Road, San-Shung
District, New Taipei City 24158 Taiwan, R.O.C.

Wellkang Ltd.,
Suite B, 29 Haley Street,
LONDON W1G 9QR,
England, United Kingdom

Importeur
SCALA Electronic GmbH
Ruhlsdorfer Straße 95
14532 Stahnsdorf, Germany
www.scala-electronic.de




Celsius
°C



Signalton



Uhrzeit



Raum-
Temperatur



Messzeit
2sec



Speicher

Ohr- und Stirn-Thermometer SC 53 FH

- ▶ Infrarot
- ▶ Messung im Ohr oder an der Stirn

4 Funktionen in einem Gerät:

- Temperaturmessung im Ohr ohne Kappe
- Temperaturmessung an der Stirn mit Kappe
- Messung der Raumtemperatur und
- Anzeige der Uhrzeit (im Ruhemodus)
- Messergebnis in ca. 2 Sek.
- Messwertspeicher (10 Werte)
- Automatischer Ruhemodus

