

PILL BOX

with vibration alarm

"Vergiss nix"



Pill Box "Vergiss nix"

- 5 medicine compartments
- set of 5 alarm times per day
- automatic daily alarm repetition
- acoustic and vibration alarm
- quartz clock with blue illuminated LCD display

The pill box "Vergiss nix" helps you to take your medication at the right time. It has five medication compartments. You can set a total of 5 different alarm times per day. These times can be set once. After midnight, the pill box automatically resets the set times for the day ahead.

As well, it is possible to set a daily recurring time for a single medicine compartment. This is useful, for example, when you only need to take birth control pills or other tablets and medications once a day.

The pill box alarm can be set in three different ways:

- acoustic alarm,
- vibration alarm,
- Acoustic AND vibration alarm.

You can easily set the alarm type using the slide switch on the side.



EAN 4022356208110

Art.-Nr.: 20811

CE
0197

More info at www.scala-electronic.de

PILL BOX

with vibration alarm

"Vergiss nix"

Specifications

Operating voltage	3V (1 x CR 2025 or CR 2032), Li content in the battery 58 mg Li, capacity 180 mAh, weight 3.0 g <u>Light and alarm functions:</u> 1,5V AAA battery (LR03), capacity 1050 mAh, weight 11 gr
Timer	Quartz timer (accuracy +/- 20 sec./month)
Alarm	optional: accoustic, vibrating or accoustic and vibrating
Alarms/day	5
Drug compartments	5
Display	LC display with 13mm characters (time)
LCD light	blue
Heart rate monitor	Number of heart beats per minute
Dimension	diameter 80 mm, height 31 mm
Device weight	96 g (incl. batteries, excl. drugs)
Colour	blue/purple
Country of origin	China
Operational manual languages	German, English, French, Dutch, Czech, Polish, Italian, Danish, Swedish, Spanish

Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.152	11.3x16.8/20.8x3.5
Sales unit	10	1.7	35x13x20
	40	7.7	43x29x37
One layer pallet	240	46.2	120x80x37
Whole pallet	1200	231	120x80x185

EAN 4022356208110

Art.-Nr.: 20811

CE
0197

More info at www.scala-electronic.de