

## PILL BOX

with vibration alarm

Vergiss nix-6



### Pill Box "Vergissnix-6"

- 6 medicine compartments
- set of 5 alarm times per day
- automatic daily alarm repetition
- acoustic and vibration alarm
- quartz clock with LCD display
- incl. timer

The pill box "Vergiss nix-6" helps you to take your medication at the right time. It has six medication compartments. You can set a total of 5 different alarm times per day. These times can be set once. After midnight, the pill box automatically resets the set times for the day ahead.



EAN 4022356208226

Art.-Nr.: 20822



0197

More info at [www.scala-electronic.de](http://www.scala-electronic.de)

## PILL BOX

with vibration alarm

Vergiss nix-6

### Specifications

Operating voltage	3V (1 x CR 2032), Li content in the battery 58 mg Li, capacity 180 mAh, weight 3.0 g
Alarm	acoustic and vibrating
Alarms/day	5
Drug compartments	6
Display	LC display
Dimension	90mmX65mmX20mm
Device weight	60g (incl. batteries, excl. drugs)
Colour	white
Country of origin	China
Operational manual languages	German, English, French, Dutch

EAN 4022356208226

Art.-Nr.: 20822



0197

More info at [www.scala-electronic.de](http://www.scala-electronic.de)