

VETERINARY THERMOMETER

for large animals

SC

1091



Veterinary thermometer SC 1091

- Temperature reading time: 15 to 25 seconds
- Measuring range: 32.0 ° C to 43.9 ° C
- Waterproof
- Long flexible measuring tip
- Big LCD-screen

The veterinary thermometer was specially developed for measuring temperature in mammals. The measuring range extends up to 43.9 ° C.

As a special feature, the thermometer has an eyelet to attach a strap. If the strap is secured around the wrist, the thermometer will not slip out of the hand even if the animal is restless or if the animal escapes.



EAN 4022356010911

Art.-Nr.: 01091

CE
0197

More info at www.scala-electronic.de

VETERINARY THERMOMETER

for large animals

Specifications



1091

Type	Maximum - Thermometer
Measuring range	32.0°C to 43.9°C Temperature below 32.0°C: display "L" for low. Temperature above 43.9°C: display "H" for high.
Measuring duration	appr. 15 to 25 seconds
Resolution	0.1°C
Deviation	± 0.1°C between 35.5°C and 42.0°C at 18°C to 28°C ambient temperature. ± 0.2°C between 32.0°C and 35.5°C and between 42.0°C and 43.9°C at 18°C to 28°C ambient temperature.
Self detection	Automatic internal check of test value 37.0°C, in case of deviation more than 0.1°C, "Err" will be displayed.
Display	LCD-Display, 3 digits
Storage temperature	Temperatures: 0°C to 50°C Relative humidity: 15% to 95%
Shutdown	appr. 3 minutes after the end of measurement
Battery info	1 x LR41 (1.5V), capacity 30mAh, weight 0.6gr
Battery indicator	yes
Water resistant	yes
Colour	white with blue details
Dimension	182.6x31.3x14.9mm
Device weight	43 g
Country of origin	China



Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0,106	23,5 x 7 x 4
Sales unit	80	8,48	41 x 60 x 27
One layer pallet	320	33,92	120 x 80 x 27
Whole pallet	1600	169,6	120 x 80 x 150

EAN 4022356010911

Art.-Nr.: 01091



0197

More info at www.scala-electronic.de

Rev.: 29.06.2022

Certified product according to the EU directive 93/42 / EEC on medical products