

SOLAR CLINICAL THERMOMETER

SC 1480



Solar Clinical Thermometer SC 1480

- Temperature reading in 10 sec,
- Rechargeable battery powered by solar energy,
- Last temperature value record
- Beep indicator
- Waterproof

The SC 1480 is a fast solar clinical thermometer, i.e. the temperature measurement is displayed after approx. 10 seconds. The final measurement is indicated by a beep signal.

With memory function, the last measured value is saved. To view the last reading, press and hold the power button for at least three seconds. After releasing the power button, the thermometer is ready for use and measure the temperature.



EAN 4022356148003

Art.-Nr.: 14800

CE
0197

more info at www.scala-electronic.de

SOLAR CLINICAL THERMOMETER

SC 1480

Specifications

Type	Solar Clinical Thermometer
Measuring range	32.0°C to 43.9°C Temperature below 32.0°C: display "L" for low. Temperature above 43.9°C: display "H" for high.
Measuring duration	appr. 10 seconds
Resolution	0.1°C
Deviation	± 0.1°C between 35.5°C and 42.0°C at 18°C to 28°C ambient temperature. ± 0.2°C between 32.0°C and 35.5°C and between 42.0°C and 43.9°C at 18°C to 28°C ambient temperature.
Fever alarm	Above 37,5°C (display is flashing)
Display	LCD-Display, 3 digits
Shutdown	appr. 3 minutes after the end of measurement
Battery info	4 solar cells / intern rechargeable battery
Battery life	approx. 40 measurements per 60 min. recharge
Battery indicator	yes
Auto detection error	display shows "ERR"
Water resistant	yes
Colour	white
Dimension	141.2 x 21.9 x 13.0 mm
Device weight	appr. 17.4 gr
Country of origin	China



	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0,025	19 x 7 x 2
Sales unit	100	6,91	28,5 x 35,5 x 42
One layer pallet	800	55,28	120 x 80 x 42
Whole pallet	3200	221,12	120 x 80 x 160

EAN 4022356148003

Art.-Nr.: 14800

CE
0197

more info at www.scala-electronic.de