

DIGITAL CLINICAL THERMOMETER

SC 1493 flex



Digital Clinical Thermometer SC 1493 flex

- flexible measuring tip
- temperature reading in 10 secs
- last temperature value record
- very sensitive unit
- compact, accurate and durable LSI (large scale integration) unit
- ABS resin body presents no hazard compared to broken mercury /glass thermometers
- fever alarm
- beep
- waterproof

The SC 1493 flex is a fast-maximum thermometer, i.e. the temperature measurement is displayed after approx. 10 seconds. The final measurement is indicated by a beep signal.

The SC 1493 flex has an integrated fever alarm. At temperatures above 37.5 ° C, the thermometer is producing short signal tones in succession.

With memory function, the last measured value is saved and auto displayed after power on. The thermometer is made by ABS resin, which is safe to use compared to broken mercury / glass thermometers.

The thermometer comes within a plastic case inside the product packaging box.



EAN 4022356014933
Art.-Nr.: 01493

CE
0197

More info at www.scala-electronic.de

Rev.: 30.11.2021

Certified product according to the EU directive 93/42 / EEC on medical products

DIGITAL CLINICAL THERMOMETER

SC 1493 flex

Specifications

Type	Maximum-Thermometer
Measuring range	32.0°C to 43.9°C Temperature below 32.0°C: display "L" for low. Temperature above 43.9°C: display "H" for high.
Measuring duration	appr. 10 seconds
Resolution	0.1°C
Deviation	± 0.1°C between 37.0°C and 39.0°C at 18°C to 28°C ambient temperature. ± 0.2°C between 35.8°C and 37.0°C and between 39.0°C and 41.0°C at 18°C to 28°C ambient temperature. ± 0.3°C for other temperature ranges.
Display	Jumbo LCD-Display, 3 digits
Shutdown	appr. 3 minutes after the end of measurement
Working temperature	Temperatures: 10°C to 40°C Relative humidity: 15% to 95%
Storage temperature	Temperatures: -25°C to 55°C Relative humidity: 15% to 95%
Battery info	1 x LR41, capacity 30mAh, weight 0.6 gr
Battery indicator	When "battery empty" symbol appears, battery must change
Colour	white with red details
Dimension	130 x 22 x 13 mm
Device weight	appr. 20 gr including battery
Country of origin	China
Operational manual languages	German, English, French, Dutch, Italian, Spanish



Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.045	14/17.5 x 5.5 x 1.9
Sales unit	10	0.56	29.5 x 11 x 6.5
	100	5.6	32 x 25x 37
One layer pallet	700	39.2	120x80x37
Whole pallet	3500	210	120x80x185

EAN 4022356014933
Art.-Nr.: 01493

CE
0197

More info at www.scala-electronic.de

Rev.: 30.11.2021

Certified product according to the EU directive 93/42 / EEC on medical products