

## DIGITAL THERMOMETER

SC 1501



### Digital Thermometer SC 1501


- flexible measuring tip
- temperature reading in 10 secs
- last temperature value record
- fever alarm
- beep
- waterproof.

The SC 1501 is a fast-maximum thermometer, i.e. the temperature measurement is displayed after approx. 10 seconds. The final measurement is indicated by a beep signal.

The SC 1501 has an integrated fever alarm. At temperatures above 37.5 ° C, the thermometer is producing short signal tones in succession.

With memory function, the last measured value is saved.

The thermometer comes within a plastic case inside the product packaging box.



EAN 4022356150105  
Art.-Nr.: 15010

CE  
0197

More info at [www.scala-electronic.de](http://www.scala-electronic.de)

# DIGITAL THERMOMETER

## Specifications

SC 1501

Measuring range	32°C ~ 43.9°C
Measuring duration	appr. 10 seconds
Resolution	0.1°C
Deviation	+/-0.1°C for temperatures 35.5°C to 42.0°C +/- 0.2°C in other areas with ambient temperature from 18°C to 28°C.
Display	LCD-Display, 3 digits
Storage temperature	Temperatures: -10°C to 50°C Relative humidity: 15% to 95%
Battery info	1 x LR41(1.5V), capacity 30mAh, weight 0.6gr
Battery indicator	When "battery empty" symbol appears, battery must change
Shutdown	Automatic after 3 min of inactivity
Dimension	130x 20 x14 mm
Device weight	appr. 13 gr including battery
Colour	white and red
Country of origin	China
Operational manual languages	German, English, French, Dutch, Czech, Polish, Italian, Danish, Swedish, Spanish, Slovakian



## Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.05	14.3/17.6x1.8x5.5 cm
Sales unit	10	0.5	30x6.5x10.5 cm
	100	6	37x23x33 cm
One layer pallet	600	36	120x80x33 cm
Whole pallet	3000	180	120x80x165 cm

EAN 4022356150105

Art.-Nr.: 15010

CE  
0197

More info at [www.scala-electronic.de](http://www.scala-electronic.de)