

DIGITAL THERMOMETER (BLUE)

SC 17



Digital thermometer SC 17 (Blue)

- temperature reading in 60 sec
- last temperature value record
- beep indicator
- waterproof

The SC 17 is a maximum thermometer, i.e. the temperature measurement is displayed after approx. 60 seconds. The final measurement is indicated by a beep signal.

With memory function, the last measured value is saved. To view the last reading, press and hold the power button for at least three seconds. After releasing the power button, the thermometer is ready for use and measure the temperature.



EAN 4022356170028

Art.-Nr.: 01048

CE
0197

more info at www.scala-electronic.de

DIGITAL THERMOMETER (BLUE)

SC 17

Specifications

Measuring range	32°C – 43.9°C
Measuring duration	appr. 60 seconds
Resolution	0.1°C
Deviation	+/- 0.1°C for temperatures 35.5°C to 42.0°C +/- 0.2°C in other areas with ambient temperature from 18°C to 28°C.
Display	LCD-Display, 3 digits
Working temperature	Temperatures: 10°C to 40°C Relative humidity: 15% to 95%
Storage temperature	Temperatures: 0°C to 50°C Relative humidity: 15% to 95%
Battery info	1 x LR41(1.5V), capacity 30mAh, weight 0.6gr
Battery indicator	When "battery empty" symbol appears, battery must change
Shutdown	Automatic after 3 min of inactivity
Dimension	126.5x 19.9 x 9.3 mm
Device weight	appr. 11.6 gr including battery
Colour	blue
Country of origin	China
Operational manual languages	German, English, French, Dutch, Czech, Polish, Italian, Danish, Swedish, Spanish, Romanian



Logistics details

	Quantity in pcs / kg	Nett weight (Kg)	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.025	0.05	17.6x5.5x1.8 cm
Sales unit	10	0.25	0.5	30x9x5.3 cm
	100	2.5	5.4	37x23x33 cm
One layer pallet	1000	25	54	120x80x33 cm
Whole pallet	5000	125	290	120x80x194 cm

EAN 4022356170028

Art.-Nr.: 01048

CE
0197

more info at www.scala-electronic.de