

DIGITAL THERMOMETER

SC 2050



Digital thermometer SC 2050

- flexible measuring tip
- colour change display (red, yellow, green)
- temperature reading in 10 secs
- last temperature value record
- fever alarm
- beep
- waterproof

The SC 2050 is a fast-maximum thermometer, i.e. the temperature measurement is displayed after approx. 10 seconds. The final measurement is indicated by a beep signal.

The SC 2050 has an integrated fever alarm. When the measured body temperature is 37.5 ° C, the fever alarm sounds (a rapid beep that sounds continuously for about 8 seconds). The LCD will flash until the thermometer is turned off.

With memory function, the last measured value is saved.

The display change colour (red, yellow or green) depending on the measured temperature.



EAN 4022356014889

Art.-Nr.: 01488

CE
0197

More info at www.scala-electronic.de

Rev.: 26.05.2021

Certified product according to the EU directive 93/42 / EEC on medical products

DIGITAL THERMOMETER

SC 2050

Specifications

Measuring range	32.0°C ~ 43.9°C
Measuring duration	appr. 10 second
Resolution	0.1°C
Deviation	+/-0.1°C for temperatures 35.5°C to 42.0°C +/- 0.2°C in other areas with ambient temperature from 18°C to 28°C.
Display	LCD-Display, 3 digits Backlight
Working temperature	Temperatures: 16°C to 35°C Relative humidity: 15% to 95%
Storage temperature	Temperatures:-10°C to 50°C Relative humidity: 15% to 95%
Battery info	1 x CR2032(3V), capacity 180mAh, weight 3gr
Battery indicator	When "battery empty" symbol appears, battery must changes
Shutdown	Automatic after 30 sec of inactivity
Dimension	117x 30 x 16 mm
Device weight	appr. 20 gr including battery
Colour	White with light blue details
Country of origin	China
Operational manual languages	German, English, French, Dutch, Czech, Polish, Italian, Danish, Swedish, Spanish



Logistics details

	Quantity in pcs / kg	Nett weight (Kg)	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.024	0.053	13.5/18x2.4x8 cm
Sales unit	10	0.24	0.6	28x9.2x13.5 cm
	100	2.4	6.7	48.5x28.5x30.3 cm
One layer pallet	400	9.6	26.8	120x80x30 cm
Whole pallet	2000	48	134	120x80x150 cm

EAN 4022356014889
Art.-Nr.: 01488

CE
0197

More info at www.scala-electronic.de

Rev.: 26.05.2021

Certified product according to the EU directive 93/42 / EEC on medical products