

VETERINARY THERMOMETER

SC

212

for large animals



Veterinary thermometer SC 212

- Temperature reading time: approx. 15 seconds
- Measuring range: 32.0 ° C to 43.9 ° C
- Waterproof
- Long measuring tip

The veterinary thermometer was specially developed for measuring temperature in mammals. The measuring range extends up to 43.9 ° C.

The special design of the thermometer enables good and quick readability even under adverse conditions.



EAN 4022356014810

Art.-Nr.: 01481

CE

0197

More info at www.scala-electronic.de

VETERINARY THERMOMETER

SC 212

for large animals

Specifications

Type	Maximum-Thermometer
Measuring range	32.0°C to 43.9°C Temperature below 32.0°C: display "L" for low. Temperature above 43.9°C: display "H" for high.
Measuring duration	appr. 15 seconds
Resolution	0.1°C
Deviation	± 0.1 ° C from 35.5 ° C to 42.0 ° C at an ambient temperature of 18 ° C to 28 ° C.
Display	LCD-Display, 3 digits
Self-test	Automatic internal check of test value 37.0°C. In case of deviation ≥ 0.1°C, "Err" will be displayed.
Shutdown	appr. 3 minutes after the end of measurement
Storage temperature	Temperatures: 0°C to 50°C Relative humidity: 15% to 95%
Battery info	1 x AAA (LR03), capacity 1050mAh, weight 11gr
Battery indicator	When "battery empty" symbol appears, battery must change
Colour	white
Device weight	appr. 56 gr including battery
Country of origin	China
Operational manual languages	German, English, French, Dutch



Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.78	6 x 20/23.5 x 2
Sales unit	12	1.008	41.5 x 13 x 7.2
	120	10.84	44 x 30 x 40
One layer pallet	720	65.04	120 x 80 x 40
Whole pallet	2880	260.16	120 x 80 x 160

EAN 4022356014810

Art.-Nr.: 01481

CE

0197

More info at www.scala-electronic.de