

DIGITAL THERMOMETER (WHITE)

SC 28 flex



Digital Thermometer SC 28 flex (white)

- flexible measuring tip
- temperature reading in 10 secs
- last temperature value record
- fever alarm
- beep
- waterproof

The SC 28 flex is a fast-maximum thermometer, i.e. the temperature measurement is displayed after approx. 10 seconds. The final measurement is indicated by a beep signal.

The SC 28 flex has an integrated fever alarm. At temperatures above 37.5 ° C, the thermometer is producing short signal tones in succession.

With memory function, the last measured value is saved and auto displayed after power on.

The thermometer comes within a plastic case inside the product packaging box.



EAN 4022356010492

Art.-Nr.: 010498



0197

More info at www.scala-electronic.de

Rev.: 25.03.2021

Certified product according to the EU directive 93/42 / EEC on medical products

DIGITAL THERMOMETER (WHITE)

SC 28 flex

Specifications

Measuring range	32°C ~ 43.9°C
Measuring duration	appr. 10 seconds
Resolution	0.1°C
Deviation	+/- 0.1°C for temperatures 35.5°C to 42.0°C +/- 0.2°C in other areas with ambient temperature from 18°C to 28°C.
Display	LCD-Display, 3 digits
Working temperature	Temperatures: 10°C to 40°C Relative humidity: 15% to 95%
Storage temperature	Temperatures: 0°C to 50°C Relative humidity: 15% to 95%
Battery info	1 x LR41(1.5V), capacity 30mAh, weight 0.6gr
Battery indicator	When "battery empty" symbol appears, battery must change
Shutdown	Automatic after 3 min of inactivity
Dimension	129x 21 x13 mm
Device weight	appr. 14 gr including battery
Colour	white
Country of origin	China
Operational manual languages	German, English, French, Dutch, Czech, Polish, Italian, Danish, Swedish, Spanish, Romanian, Croatian, Russian



Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.05	14.2x5.5x1.8 cm
Sales unit	14	0.813	29x15x7 cm
	140	8.64	37.5x32.5x31 cm
One layer pallet	720	68.8	120x80x45 cm
Whole pallet	4320	287.8	120x80x175 cm

EAN 4022356010492

Art.-Nr.: 010498



0197

More info at www.scala-electronic.de

Rev.: 25.03.2021

Certified product according to the EU directive 93/42 / EEC on medical products