

VETERINARY THERMOMETER

for small animals

SC

312



Veterinary thermometer SC 312

- Temperature reading time: approx. 20 seconds
- Measuring range: 32.0 ° C to 43.9 ° C
- Waterproof
- Standard measuring tip

The veterinary thermometer was specially developed for temperature measurement in mammals. The measuring range extends up to 43.9 ° C.

As a special feature, the thermometer has an eyelet to attach a strap. If the strap is secured around the wrist, the thermometer will not slip out of the hand even if the animal is restless or if the animal escapes.



EAN 4022356014827

Art.-Nr.: 01482

CE
0197

More info at www.scala-electronic.de

VETERINARY THERMOMETER

for small animals

Specifications

Type	Maximum-Thermometer
Measuring range	32.0°C to 43.9°C Temperature below 32.0°C: display "L" for low. Temperature above 43.9°C: display "H" for high.
Measuring duration	appr. 20 seconds
Resolution	0.1°C
Deviation	± 0.1°C between 35.0°C and 42.0°C at 18°C to 28°C ambient temperature ± 0.2°C for other temperature ranges at 18°C to 28°C ambient temperature
Display	LCD-Display, 3 digits
Self-test	Automatic internal check of test value 37.0°C. In case of deviation ≥ 0.1°C, "Err" will be displayed.
Shutdown	appr. 3 minutes after the end of measurement
Working temperature	Temperatures: 10°C to 40°C Relative humidity: 15% to 95%
Storage temperature	Temperatures: 0°C to 50°C Relative humidity: 15% to 95%
Battery info	1 x LR41, capacity 30mAh, weight 0.6 gr
Battery indicator	When "battery empty" symbol appears, battery must change
Colour	white with blue details
Dimension	124.4x 18.3 x 9.7mm
Device weight	appr. 11.4 gr including battery
Country of origin	China
Operational manual languages	German, English, French, Italian

SC 312



Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.042	7x14.8/18.2x1.5
Sales unit	12	0.557	30x8x11
	120	6.15	45x33x24
One layer pallet	960	49.2	120x80x24
Whole pallet	6720	344.4	120x80x170

EAN 4022356014827

Art.-Nr.: 01482

CE
0197

More info at www.scala-electronic.de