

## DIGITAL BODY THERMOMETER - RED

SC 37T



### Digital Body Thermometer SC 37T - Red

- temperature reading in 60 sec
- last temperature value record
- beep indicator
- waterproof

The SC 37T is a maximum thermometer, i.e. the temperature measurement is displayed after approx. 60 seconds. The final measurement is indicated by a beep signal.

With memory function, the last measured value is saved. To view the last reading, press and hold the power button for at least three seconds. After releasing the power button, the thermometer is ready for use and measure the temperature.

The thermometer comes within a plastic case inside the product packaging box.



EAN 4022356014414

Art.-Nr.:01441

CE  
0197

More info at [www.scala-electronic.de](http://www.scala-electronic.de)

Rev.: 02.06.2021

Certified product according to the EU directive 93/42 / EEC on medical products

## DIGITAL BODY THERMOMETER - RED

SC 37T

### Specifications

Measuring range	32 <sup>0</sup> C ~ 43.9 <sup>0</sup> C
Measuring duration	appr. 60 seconds
Resolution	0.1 <sup>0</sup> C
Deviation	+/- 0.1 <sup>0</sup> C for temperatures 35.5 <sup>0</sup> C to 42.0 <sup>0</sup> C +/- 0.2 <sup>0</sup> C in other areas with ambient temperature from 18 <sup>0</sup> C to 28 <sup>0</sup> C.
Display	LCD-Display, 3 digits
Working temperature	Temperatures: 10 <sup>0</sup> C to 40 <sup>0</sup> C Relative humidity: 15% to 95%
Storage temperature	Temperatures: 0 <sup>0</sup> C to 50 <sup>0</sup> C Relative humidity: 15% to 95%
Battery info	1 x LR41(1.5V), capacity 30mAh, weight 0.6gr
Battery indicator	When "battery empty" symbol appears, battery must change
Shutdown	Automatic after 3 min of inactivity
Dimension	124.4x 18.3 x 9.7 mm
Device weight	appr. 11.4 gr including battery
Colour	Red
Country of origin	China
Operational manual languages	German, English, French, Dutch, Czech, Polish, Italian, Danish, Swedish, Spanish, Romanian, Turkish, Russian



### Logistics details

	Quantity in pcs / kg	Nett weight (Kg)	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.025	0.045	17.8x1.2x7.2 cm
Sales unit	12	0.3	0.6	21.5x7.5x19 cm
	120	3	6.4	40x23.5x39 cm
One layer pallet	1080	27	57.6	120x80x39 cm
Whole pallet	5400	135	288	120x80x195 cm

EAN 4022356014414

Art.-Nr.:01441

CE  
0197

More info at [www.scala-electronic.de](http://www.scala-electronic.de)

Rev.: 02.06.2021

Certified product according to the EU directive 93/42 / EEC on medical products