

DIGITAL MEDICAL THERMOMETER

Display sorted by colour (24 pieces)

SC 42TM



Digital Thermometer SC 42TM



- flexible measuring tip
- temperature reading in 10 secs
- last temperature value record
- fever alarm
- beep
- waterproof

The SC 42TM is a fast-maximum thermometer, i.e. the temperature measurement is displayed after approx. 10 seconds. The final measurement is indicated by a beep signal.

The SC 42TM has an integrated fever alarm. At temperatures above 37.5 °C, the thermometer is producing short signal tones in succession.

With memory function, the last measured value is saved and auto displayed after power on.

The thermometer comes within a plastic case inside the product packaging box. There are 2 colours available (red and blue) to choose from - 12 pieces per colour sorted in the display packaging box.



EAN- 4022356014490

Art.-Nr.: 01452

CE
0197

More info at www.scala-electronic.de

Rev.:25.03.2021

Certified product according to the EU directive 93/42 / EEC on medical products

DIGITAL MEDICAL THERMOMETER

Display sorted by colour (24 pieces)

SC 42TM

Specifications

Measuring range	32°C ~ 43.9°C
Measuring duration	appr. 10 seconds
Resolution	0.1°C
Deviation	+/- 0.1°C for temperatures 35.5°C to 42.0°C +/- 0.2°C in other areas with ambient temperature from 18°C to 28°C.
Display	LCD-Display, 3 digits
Working temperature	Temperatures: 10°C to 40°C Relative humidity: 15% to 95%
Storage temperature	Temperatures: 0°C to 50°C Relative humidity: 15% to 95%
Battery info	1 x LR41(1.5V), capacity 30mAh, weight 0.6gr
Battery indicator	When "battery empty" symbol appears, battery must changes
Shutdown	Automatic after 3 min of inactivity
Dimension	129.5x 21 x 12 mm
Device weight	appr. 12 gr including battery
Colour	White (red or blue button)
Country of origin	China
Operational manual languages	German, English, French, Dutch, Czech, Polish, Italian, Danish, Swedish, Spanish, Romanian



Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.06	20/17.5x2.2x6.3 cm
Sales unit	12	0.738	29x7x21 cm
	120	8.11	59.5x24x36 cm
One layer pallet	960	70	120x80x36 cm
Whole pallet	4800	349.4	120x80x160 cm

EAN- 4022356014490

Art.-Nr.: 01452

CE
0197

More info at www.scala-electronic.de

Rev.:25.03.2021

Certified product according to the EU directive 93/42 / EEC on medical products