

WRIST – BLOOD PRESSURE MONITOR

SC 6400



Wrist blood pressure monitor SC 6400 Blue

- WHO indicator
- beep signal during the measurement process
- memory for 85 measurements
- display of date and time

The blood pressure is measured on the wrist. The blood pressure monitor has a WHO indicator. The measured values of blood pressure are according to the guidelines of the World Health Organization.

After being switched on, the SC 6400 automatically inflates until the factory-set inflation pressure is reached. If this is not sufficient, the device detects this and automatically pumps up (Fuzzy Logic).

The subsequent measurement process is accompanied by a beep until the blood pressure values are determined and stops when the measurement is completed.

After a successful measurement, the systolic and diastolic blood pressure and the pulse are displayed at the same time.



EAN 4022356025021

Art.-Nr.: 2184

CE
0197

More info at www.scala-electronic.de

WRIST – BLOOD PRESSURE MONITOR

Specifications

Measuring range	Pressure: 20 ~ 300mmHg Pulse rate: 40 ~ 200 pulse / min
Resolution	1 mmHg
Accuracy	Pressure: ± 3 mmHg Pulse frequency: $\pm 5\%$ of the measured value
Measurement method	Oscillometric
Controls	3 buttons: On / Off (O / I), memory (M) & User button
Display	LCD display, WHO bar graph
Time	Display of date & time
Memory	85 sets for systole / diastole / pulse, date/time
Working temperature	Temperatures: 5°C to 40°C Relative humidity: 15% to 90%
Storage temperature	Temperatures: -10°C to 55°C Relative humidity: up 93%
Voltage operation	2xAAA (LR03) alkaline batteries, capacity 2x1050 mAh, Battery weight 2x11 g
Battery indicator	If the symbol for dead batteries appears on the display, replace all batteries with new ones at the same time.
Shutdown	Automatic after appr. 150 sec of inactivity
Dimension	77.8x72.3x30.7mm (without cuff)
Cuff size	for wrist circumferences from 135 to 195 mm
Device weight	appr. 108.2 gr (without batteries)
Colour	White with blue details
Country of origin	China
Operational manual languages	German, English, French, Dutch, Czech, Polish, Italian, Danish, Romanian, Spanish, Turkish, Portuguese

SC 6400



Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.235	9.5x18x10 cm
Sales unit	18	4.65	26x33x20 cm
	36	10	41x28x35 cm
One layer pallet	288	80	120x80x35 cm
Whole pallet	1440	400	120x80x175 cm

EAN 4022356025021

Art.-Nr.: 2184

CE
0197

More info at www.scala-electronic.de