

UPPER ARM BLOOD PRESSURE MONITOR

SC 7670



Upper arm blood pressure monitor SC 7670

- WHO indicator
- memory for 2 x 50 measurements
- beep signal during the measurement process
- display of date and time

The blood pressure is measured on the upper arm. The blood pressure monitor has a WHO indicator. The measured values of blood pressure are according to the guidelines of the World Health Organization.

After being switched on, the device automatically inflates until the factory-set inflation pressure is reached. If this is not sufficient, the device detects this and automatically pumps up (Fuzzy Logic).

The subsequent measurement process is accompanied by a beep until the blood pressure values are determined and stops when the measurement is completed.

After a successful measurement, the systolic and diastolic blood pressure and the pulse are displayed at the same time.



EAN 4022356024956

Art.-Nr.: 02495

CE
0197

more info at www.scala-electronic.de

UPPER ARM BLOOD PRESSURE MONITOR

Specifications

Measuring range	Pressure: 20 ~ 300mmHg Pulse rate: 40 ~ 200 pulse / min
Accuracy	Pressure: 3 mmHg Pulse frequency: $\pm 5\%$ of the measured value
Measurement method	Oscillometric
Controls	3 buttons: On / Off (O / I), memory (M) & user button
Display	LCD display, digital liquid crystal display, WHO bar graph
Memory	2x50 sets
Time	Display date & time
Working temperature	Temperatures: 5°C to 40°C Relative humidity: 15% to 95%
Storage temperature	Temperatures: -20°C to 55°C Relative humidity: up 93%
Battery info	4 x AA(LR6), capacity 4x2300mAh, weight 4x23 gr
Battery indicator	When "battery empty" symbol appears, battery must changes
Shutdown	Automatic after appr. 150 sec of inactivity
Dimension	142 x 142 x 120 mm
Cuff size	For upper arm circumferences from 20 to 35 cm
Device weight	appr. 522 gr including battery
Colour	White with blue display
Country of origin	China
Operational manual languages	German, English, French, Dutch, Italian



Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0.87	17.5x14.5x17.5 cm
Sales unit	4	3.98	60x19x19.5 cm
	16	15.9	61.5x40x40 cm
One layer pallet	64	63.5	120x80x40 cm
Whole pallet	256	254	120x80x160 cm

EAN 4022356024956

Art.-Nr.: 02495



more info at www.scala-electronic.de

Rev.: 28.09.2021

Certified product according to the EU directive 93/42 / EEC on medical products