

## INFRARED FOREHEAD THERMOMETER

SC 8360



### Infrared forehead thermometer SC 8360

- measurement over the forehead - contactless, or object
- 1 approx. second measuring time
- 9 memory sets
- beep
- fever alarm
- temperature display with backlight
- battery operation

With the infrared forehead thermometer SC 8360, the human body temperature can be determined quickly and easily. The measuring point is the forehead. The body temperature is determined by evaluating the infrared radiation. The device can be used on both adult patients and children.

The SC 8360 has an integrated fever alarm. The measurement takes 1 second. After a successful measurement, the last nine measured values will be saved. Batteries are included and have a life time of around 3,000 measurements.



EAN 4022356083601  
Art.-Nr.:08360

CE  
0197

# INFRARED FOREHEAD THERMOMETER

SC 8360

## Specifications

Measuring type	Infrared
Measuring location	Forehead (contactless) / Object
Measuring distance	2 cm
Measuring range	Forehead: 32.0°C - 43.3°C Object: 0 - 100°C
Measuring duration	appr. 1 second
Resolution	0.1°C
Deviation	Forehead: 35.5°C - 42.0°C: +/- 0.2°C other range +/- 0.3°C Object: +/- 1.0°C
Fever alarm	From 37.5°C
Display	LCD-Display, 3 digits, Backlight
Operational temperature	Temperatures: 16°C to 35°C, R.H. of 15% to 90% non-condensing
Storage temperature	Temperatures: -20°C to +50°C, R.H. of 15% to 90% non-condensing
Battery info	2 x (R03, LR03 (1.5V), Capacity 1050mAh, weight 11 gr)
Battery indicator	When "battery empty" symbol appears, battery must change
Memory	9 sets
Shutdown	automatically after appr. 10sec of inactivity
Dimension	41.4 x 44 x 124 mm
Device weight	appr. 88 gr including battery
Colour	white (blue button)
Certification	ISO 13485 and compliance with ISO 80601-2-56
Country of origin	China
Operational manual languages	German, English, French, Dutch, Czech, Polish, Italian, Danish, Swedish, Spanish



## Logistics details

	Quantity in pcs / kg	Gross weight (Kg)	Dimensions (cm)
Single package	1	0,137	8x16,2/20x6
Sales unit	10	1,49	34x9x33
	100	16,98	67x51x35
One layer pallet	400	67,92	120x80x35
Whole pallet	2000	339,6	120x80x175

EAN 4022356083601  
Art.-Nr.: 08360

CE  
0197

More info at [www.scala-electronic.de](http://www.scala-electronic.de)

Rev.: 31.08.2021

Certified product according to the EU directive 93/42 / EEC on medical products